

Poster 77**CEREBRAL OXIMETRY VIA NEAR INFRARED SPECTROSCOPY MIRRORS DESATURATIONS OF BRAIN TISSUE PARTIAL PRESSURE OF OXYGEN: A CASE STUDY**

Elissa Wible, Carmelo Graffagnino, DaiWai Olson, Michael James

Duke University, Durham, NC, United States

Introduction:

Near infrared spectroscopy is a non-invasive method of monitoring cerebral oxygenation. By employing time and spatial resolution of several light wavelengths, cortical blood flow, volume and oxygenation can be quantified (cerebral oximetry). We present a case utilizing cerebral oximetry in a patient with cerebral vasospasm after subarachnoid hemorrhage (SAH) with concurrent use of brain tissue partial pressure of oxygen (PbtO₂) monitoring.

Methods:

A 46 year old woman developed severe diffuse vasospasm following SAH. We monitored intracranial pressure (ICP) as well as tissue oximetry (PbtO₂) via Licox (Integra) catheter placed in the distribution of the left MCA. Over a 24-hour period during the third day of vasospasm, CerOx (Ornim) monitoring was applied over the left fronto-temporal area to evaluate the relationship between PbtO₂ and non-invasive cerebral oximetry.

Results:

There were 12 episodes of PbtO₂ desaturation (<20 mmHg for >5 mins) over the period of dual monitoring. Over 80% of these PbtO₂ desaturations were preceded by >20% decline in CerOx values from baseline. There were 22 episodes of cerebral oximetry desaturation (<50% for >10 mins). Less than 50% of CerOx desaturations were temporally related to a decline in PbtO₂ to less than 20 mmHg. Hemoglobin was stable at 10 mg/dL and ICP was well controlled (<15 mmHg) during the entire 24 hours.

Conclusions:

In this subject, desaturations of PbtO₂ appeared to be related to desaturations by non-invasive cerebral oximetry; the converse was not the case. Perhaps cortical oxygen desaturations (CerOx) occur with increased frequency compared to subcortical oxygen desaturations (PbtO₂) in diffuse vasospasm after SAH. CerOx monitoring may provide an enhanced understanding of oxygen delivery and utilization during periods of ongoing cerebral ischemia. Further studies are required to substantiate these findings.

Financial Support: None